

Pendants

H-Beam Pendants for High Care

Extruded aluminium H-Beams can be customized to any length. Used to serve multiple beds in a high-care environment. Unobstructed access from behind bed. Vertical columns can slide and twist.

Specifications:

Electrical supply: services through ceiling void/core drill hole from plant room.

- High voltage: 230 V
- Low voltage: to suit specific requirements – data, telephone, network plugs etc.
- Medical gas: the unit can accommodate any number of gas sockets.

Features:

- Slide and twist function
- Customizable to client specification: colors, gas outlets, electricals

