

Pendants

Various Theatre Pendants for Stacks and Screens

An integrated ceiling-mounted pendant, with articulated arms and body constructed from extruded aluminium, body is covered with a hard-wearing powder-coated finish. Offers a flexible “off the floor” solution that allows doctors and nurses to access the necessary services without the risk of exposed piping and electrical cables on the floor. Allows for uninterrupted maintenance

Specifications:

Electrical supply: services through ceiling void/core drill hole from plant room.

- High voltage: 230 V
- Low voltage: to suit specific requirements – data, telephone, network plugs etc.
- Medical gas: the unit can accommodate a total of eight gas sockets.
- Product weight: 87 kg
- Rotation: 320 °

Features:

- Rotation of 320 °
- Customizable to client specification: colors, gas outlets, electricals

